



APEXTM

KERATOCONUS CONTACT LENSWITH FITTING FORGIVENESS

The Apex lens design for keratoconus incorporates a number of fitting forgiveness factors to simplify fitting the keratoconus patient. These inherent design properties are intended to aid the practitioner in managing the various topographies seen with this corneal disorder.

This new proprietary lens design is based upon the belief that those patients who are successful with gas permeable lens designs must be fitted with a lens that minimizes apical interaction, and avoids impingement of the peripheral cornea, specifically in the superior quadrant.

LENS PARAMETER AVAILABILITY

Base Curve	68.00D (4.96) to 40.00D (8.44)
Diameter	7.9 to 10.5
Power	+30.00D to -30.00D
Periphery	Steep, Medium, Flat

DIAGNOSTIC SET PARAMETERS

Standart 24-Lens Set

Base Curve	58.00D (5.82) to 40.00D (8.44)	
Diameter	8.2 to 10.0	
Power	Flat to -17.00	
Periphery	Medium	
Material	Available in any material	

A wratten filter is included.

Additional 10-Lens Diagnostic Set Available

Base Curve	68.00D (4.96) to 59.00D (5.72)
Diameter	8.2
Power	-17.00 to -20.00
Periphery	Medium
Material	Available in any material





TROUBLESHOOTING

Patient Symptoms	Possible Causes	Plan
Long Diding High		Steepen Base Curve by 1.00 D
Lens Riding High		Reduce Diameter by .3mm
Lana Didina Law	Movement with Blink	Flatten Base Curve by 1.00 D
Lens Riding Low	No Movement	Decrease Diameter by .3mm
F M		Steepen Base Curve by 1.00 D
Excessive Movement		Increase Diameter by .3mm
Restricted Movement		Flatten Base Curve by 1.00 D
Nestricted Provenient		Decrease Diameter by .3mm
	Centrally	Flatten Base Curve by 1.00 D
Bubbles		Decrease Diameter by .3mm
Dubbles	Peripherally	Decrease Edge Lift (Flat to Medium, Medium to Steep)
		Increase Diameter by .3mm
	New Lens	Clean Lens
Surface Non		Change Material
Wetting	Older Lens	Clean Lens
		Order New Lens
Too Much Edge Lift		Reduce Edge Lift (Flat to Medium, Medium to Steep)
		Increase Diameter by .3 mm
Not Enough Edge Lift		Increase Edge Lift (Steep to Medium, Medium to Flat)
		Decrease Diameter by .3 mm

Apex™ Keratoconus

