

# 3-MONTHLY REPLACEMENT

## QUATTRO

### HYDROGEL



**Quattro** provides spherical, toric and multifocal correction in multiple diameters for patients already accustomed to a 3-monthly lens replacement.






### PARAMETERS

<b>Base curves (mm)</b>	SPH, TOR 7.70 to 9.80 (0.30) (Ø14.50) MF 8.00 to 9.00 (0.20) (Ø14.00) SPH, TOR 7.10 to 9.20 (0.30) (Ø13.00)
<b>Diameters (mm)</b>	SPH, TORIC 13.00 & 14.50 MF 14.00
<b>Spheres (D)</b>	SPH, TOR: ±30.00 (0.25) MF: -12.00 to -1.00 / +1.00 to +8.00 (0.25)
<b>Cylinders (D)</b>	-0.75 to -8.00 (0.25)
<b>Axes (°)</b>	All (5°)

	SPH +	SPH -
A	1.00 CN	1.00 CD
B	1.75 CN	2.00 CD
C	2.50 CN	3.00 CD

### OPTICAL DESIGN



Lens design parameters may change depending on the power

### MATERIAL

<b>Type</b>	Filcon 1 (15) [49%]
<b>Dk (iso 9913-1-1998)</b>	15
<b>DK/T (-3.00D)</b>	17
<b>Water content</b>	49%
<b>Cof</b>	0.09
<b>Modulus</b>	0.41
<b>Handling tint</b>	BLUE
<b>Pack size</b>	Single and 2-pack
<b>Manufacturing process</b>	LATHED