



**GAS PERM LENS SYSTEM WITH  
ECCENTRICITY CONTROL**

**Advanced system for keratoconus**

**#specialize #lamspectrum**


**S P E C T R U M**  
INTERNATIONAL




[www.spctinternational.com](http://www.spctinternational.com)  
[www.spectrumcuniversity.com](http://www.spectrumcuniversity.com)

Made by







Elipsys STD #1	Elipsys STD #2	Elipsys STD #3	Elipsys STD #4	Elipsys STD #5	Elipsys STD #6
BC: 8.40 Pwr: +1.50 Dia: 9.80 Ecc: STD	BC: 8.20 Pwr: +0.50 Dia: 9.80 Ecc: STD	BC: 8.00 Pwr: -0.50 Dia: 9.80 Ecc: STD	BC: 7.90 Pwr: -1.00 Dia: 9.80 Ecc: STD	BC: 7.80 Pwr: -1.50 Dia: 9.80 Ecc: STD	BC: 7.70 Pwr: -2.00 Dia: 9.80 Ecc: STD




Elipsys STD #7	Elipsys STD #8	Elipsys STD #9	Elipsys STD #10	Elipsys KC #11	Elipsys KC #12
BC: 7.60 Pwr: -2.50 Dia: 9.80 Ecc: STD	BC: 7.50 Pwr: -3.00 Dia: 9.80 Ecc: STD	BC: 7.40 Pwr: -3.50 Dia: 9.80 Ecc: STD	BC: 7.30 Pwr: -4.00 Dia: 9.80 Ecc: STD	BC: 7.60 Pwr: -2.50 Dia: 9.60 Ecc: MED	BC: 7.40 Pwr: -3.75 Dia: 9.60 Ecc: MED 2



Elipsys KC #13	Elipsys KC #14	Elipsys KC #15	Elipsys KC #16	Elipsys KC #17	Elipsys KC #18
BC: 7.20 Pwr: -5.00 Dia: 9.60 Ecc: MED 2	BC: 7.00 Pwr: -6.25 Dia: 9.60 Ecc: FLAT	BC: 6.80 Pwr: -7.75 Dia: 9.60 Ecc: FLAT	BC: 6.60 Pwr: -9.25 Dia: 9.60 Ecc: FLAT	BC: 6.40 Pwr: -10.75 Dia: 9.60 Ecc: FLAT	BC: 6.20 Pwr: -12.50 Dia: 9.60 Ecc: FLAT





Elipsys KC #19	Elipsys KC #20	Elipsys SE #21	Elipsys SE #22	Elipsys SE #23	Elipsys SE #24
BC: 6.00 Pwr: -14.25 Dia: 9.60 Ecc: FLAT 2	BC: 5.80 Pwr: -16.25 Dia: 9.60 Ecc: FLAT 2	BC: 8.50 Pwr: +2.00 Dia: 11.20 Ecc: STD	BC: 8.30 Pwr: +1.00 Dia: 11.20 Ecc: STD	BC: 8.10 Pwr: +0.00 Dia: 11.20 Ecc: STD	BC: 7.90 Pwr: -1.00 Dia: 11.20 Ecc: MED




Elipsys SE #25	Elipsys SE #26	Elipsys SE #27	Elipsys SE #28	Elipsys SE #29	Elipsys SE #30
BC: 7.70 Pwr: -2.00 Dia: 11.20 Ecc: MED	BC: 7.50 Pwr: -3.25 Dia: 11.20 Ecc: MED 2	BC: 7.30 Pwr: -4.50 Dia: 11.20 Ecc: FLAT	BC: 7.10 Pwr: -5.75 Dia: 11.20 Ecc: FLAT	BC: 6.90 Pwr: -7.25 Dia: 11.20 Ecc: FLAT	BC: 6.70 Pwr: -8.75 Dia: 11.20 Ecc: FLAT 2

Aligned snap closure to secure closure airtight by means of a security strip.

Case attached to the base	Remove security strip
 <p>A presión</p>	

**TEST AREA**



**OD**                      **OS**