



SISTEMA DE LENTE DE CONTACTO GP
CON CONTROL DE EXCENTRICIDAD Y
APOYO PERIFÉRICO

Sistema avanzado para Queratocono

#Especialízate #YoSoyspectrum

S P E C T R U M
I N T E R N A T I O N A L

www.spcinternational.com

www.spectrumcuniversity.com




[/spctinternational](https://www.facebook.com/spctinternational)

Hecho por


Lenticon
EUROPE




• PREMIUM QUALITY •




Elipsys STD #1	Elipsys STD #2	Elipsys STD #3	Elipsys STD #4	Elipsys STD #5	Elipsys STD #6
BC: 8.40 Pwr: +1.50 Dia: 9.80 Exc: STD	BC: 8.20 Pwr: +0.50 Dia: 9.80 Exc: STD	BC: 8.00 Pwr: -0.50 Dia: 9.80 Exc: STD	BC: 7.90 Pwr: -1.00 Dia: 9.80 Exc: STD	BC: 7.80 Pwr: -1.50 Dia: 9.80 Exc: STD	BC: 7.70 Pwr: -2.00 Dia: 9.80 Exc: STD




Elipsys STD #7	Elipsys STD #8	Elipsys STD #9	Elipsys STD #10	Elipsys KC #11	Elipsys KC #12
BC: 7.60 Pwr: -2.50 Dia: 9.80 Exc: STD	BC: 7.50 Pwr: -3.00 Dia: 9.80 Exc: STD	BC: 7.40 Pwr: -3.50 Dia: 9.80 Exc: STD	BC: 7.30 Pwr: -4.00 Dia: 9.80 Exc: STD	BC: 7.60 Pwr: -2.50 Dia: 9.60 Exc: MED	BC: 7.40 Pwr: -3.75 Dia: 9.60 Exc: MED 2



Elipsys KC #13	Elipsys KC #14	Elipsys KC #15	Elipsys KC #16	Elipsys KC #17	Elipsys KC #18
BC: 7.20 Pwr: -5.00 Dia: 9.60 Exc: MED 2	BC: 7.00 Pwr: -6.25 Dia: 9.60 Exc: FLAT	BC: 6.80 Pwr: -7.75 Dia: 9.60 Exc: FLAT	BC: 6.60 Pwr: -9.25 Dia: 9.60 Exc: FLAT	BC: 6.40 Pwr: -10.75 Dia: 9.60 Exc: FLAT	BC: 6.20 Pwr: -12.50 Dia: 9.60 Exc: FLAT



Elipsys KC #19	Elipsys KC #20	Elipsys SE #21	Elipsys SE #22	Elipsys SE #23	Elipsys SE #24
BC: 6.00 Pwr: -14.25 Dia: 9.60 Exc: FLAT 2	BC: 5.80 Pwr: -16.25 Dia: 9.60 Exc: FLAT 2	BC: 8.50 Pwr: +2.00 Dia: 11.20 Exc: STD	BC: 8.30 Pwr: +1.00 Dia: 11.20 Exc: STD	BC: 8.10 Pwr: +0.00 Dia: 11.20 Exc: STD	BC: 7.90 Pwr: -1.00 Dia: 11.20 Exc: MED




Elipsys SE #25	Elipsys SE #26	Elipsys SE #27	Elipsys SE #28	Elipsys SE #29	Elipsys SE #30
BC: 7.70 Pwr: -2.00 Dia: 11.20 Exc: MED	BC: 7.50 Pwr: -3.25 Dia: 11.20 Exc: MED 2	BC: 7.30 Pwr: -4.50 Dia: 11.20 Exc: FLAT	BC: 7.10 Pwr: -5.75 Dia: 11.20 Exc: FLAT	BC: 6.90 Pwr: -7.25 Dia: 11.20 Exc: FLAT	BC: 6.70 Pwr: -8.75 Dia: 11.20 Exc: FLAT 2

Cierre a presión alineado para asegurar el cierre hermético por medio de una tira de seguridad.

Estuche acoplado a la base	Remover tira de seguridad
	

ÁREA DE PRUEBAS



OD **OS**