



The SPECTRUM SALINE Isotonic buffered saline **solution** is a sterile, isotonic and buffered solution in single-dose format without preservatives that acts by mechanical entrainment, which is why it is especially suitable for people with sensitive eyes and nostrils.

FREE PRESERVATIVE

Which are your principal functions?

With contact lenses:



- Rinse and insert scleral lenses.
- Rinse, hydrate and preserve scleral, soft and rigid gas permeable lenses.



Ocular instillation:

It is a very useful complementary solution for users of all types of contact lenses since it eliminates, for example, possible dust particles, reduces the discomfort derived from adverse environmental conditions, such as environments with air conditioning / heating, dry environments, environments with smoke., etc.

Technical data:

SPECTRUM SALINE Isotonic buffered saline solution is composed of elements that make it very effective:

- Sodium Chloride: For tonicity and osmolarity.
- Boric Acid/Sodium Borate Buffer Control the pH of the solution to be the pH of the eye.
- Single-dose format without preservatives

